



Recombinant Chicken Argininosuccinate synthase (ASS1)

Product Code	CSB-MP730490CH
Abbreviation	ASS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZJ23
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAEPRDTVVL AYSGGLDTSC ILVWLKEQGY TVIAFLANIG QTEDFDAAQK KALALGAKKV YIQDVCREFV EDFIWPAVRA NALYEDRYML GSALARPCIA RHLVLIAQEE GARYIAHGAT GKGNDQVRFE LGCYALCPSI KVIAPWRMPE FYQRFPGRR LMEYAQKHGI PVPVTPKAPW SMDENLMHIS YEAGILENPK NRAPLDLYTK TCNPTTSPDV PDELEIEFEK GVPVKVTNTR NGVTHRSAL LFVYLNDIAS KHGVGRVDIV ENRFVGMKSR GIYETPAGTI LYHAHLIDIEA FTMDREVRKI KQGLSLKFSE LVYNGFWYSP ECEFLKHICIA RSQQAVAGTV HVSVFKGQVY VLGRESPHSL YNEELVSM DV QGDYEPADAT GFININALRL KEYHRLQSKV STKQDE
Source	Mammalian cell
Target Names	ASS1
Protein Names	Recommended name: Argininosuccinate synthase EC= 6.3.4.5 Alternative name(s): Citrulline--aspartate ligase
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein catalyzes the penultimate step of the arginine biosynthetic pathway. There are approximately 10 to 14 copies of this gene including the pseudogenes scattered across the human genome, among which the one located on chromosome 9 appears to be the only functional gene for argininosuccinate synthetase. Mutations in the chromosome 9 copy of ASS cause citrullinemia. Two transcript variants encoding the same protein have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.