



# Recombinant Chicken Argininosuccinate synthase (ASS1)

<b>Product Code</b>	CSB-YP730490CH
<b>Abbreviation</b>	ASS1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5ZJ23
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAEPRDTVVL AYSGGLDTSC ILVWLKEQGY TVIAFLANIG QTEDFDAAQK KALALGAKKV YIQDVCREFV EDFIWPAVRA NALYEDRYML GSALARPCIA RHLVLIAQEE GARYIAHGAT GKGNDQVRFE LGCYALCPSI KVIAPWRMPE FYQRFPGRRE LMEYAQKHGI PVPVTPKAPW SMDENLMHIS YEAGILENPK NRAPLDLYTK TCNPTTSPDV PDELEIEFEK GVPVKVTNTR NGVTHRSALE LFVYLNDIAS KHGVGRVDIV ENRFVGMKSR GIYETPAGTI LYHAHLDIEA FTMDREVRKI KQGLSLKFSE LVYNGFWYSP ECEFLKHICIA RSQQAVAGTV HVSVFKGQVY VLGRESPHSL YNEELVSM DV QGDYEPADAT GFININALRL KEYHRLQSKV STKQDE
<b>Source</b>	Yeast
<b>Target Names</b>	ASS1
<b>Protein Names</b>	Recommended name: Argininosuccinate synthase EC= 6.3.4.5 Alternative name(s): Citrulline--aspartate ligase
<b>Expression Region</b>	1-416
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein catalyzes the penultimate step of the arginine biosynthetic pathway. There are approximately 10 to 14 copies of this gene including the pseudogenes scattered across the human genome, among which the one located on chromosome 9 appears to be the only functional gene for argininosuccinate synthetase. Mutations in the chromosome 9 copy of ASS cause citrullinemia. Two transcript variants encoding the same protein have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.