



Recombinant Chicken WD repeat-containing protein 61 (WDR61)

Product Code	CSB-MP730495CH
Abbreviation	WDR61
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZJH5
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MTTQYGILFK QEQAHDDAIW SVAWGKNKKN GSETVISGSL DDLVKVWKWN DEKLDLQWTL EGHQLGVVSV DISHTGSIAA SSSLDAHIRL WDLETGKQIK SIDAGPVD AW SLAFSPDSQF LATGSHVGKV NIFGVETGKK EYSLDTRGKF ILSIAYSPDG KYLASGAIDG IINIFDIATG KLLHTLEGHA MPIRSLTFSP DSQLLV T ASD DGYIKIYDVQ HANLAGT LSG HGSWVLN VAF CPDDTHFVSS SSDKSVK VWD AGTRTCVHTF FDHQDQVWGV KYNGSGSKIV SVGDDQEIH YDCPV
Source	Mammalian cell
Target Names	WDR61
Protein Names	Recommended name: WD repeat-containing protein 61
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.