



Recombinant Chicken Protein FAM60A (FAM60A)

Product Code	CSB-YP730502CH
Abbreviation	FAM60A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZJV7
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MFGFRKPKMY RSIEGCCICR AKSSSSRFTD SKRYEKDFQN CFGLHEARSG DICNACVLLV KRWKKLPAGS KKNWNHVVDA RAGPSLKTTL KPKKMKTLTG SRIKSNQISK LQKEFKRHNS DAHSTTSSAS PAQSPCYSNQ SDEGSDTEMS AGSSRTPVFS FLDLTYWKRQ KVCCGIIYKG RFGEVLIDTH LFKPCCSNKK SATEKPEQEG PQSPAISTQE EW
Source	Yeast
Target Names	SINHCAF
Protein Names	Recommended name: Protein FAM60A
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.