



Recombinant Legionella pneumophila subsp. pneumophila Probable rRNA maturation factor (lpg1440)

Product Code	CSB-MP730554LDL
Abbreviation	lpg1440
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZVK0
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)
Purity	≥85% (SDS-PAGE)
Sequence	MTYHIDIQNA TGKLLPLSED EITKLASLAL RDHKQDAELT VRLVDVEEMT YLNHTYRKKN KPTNVLAFCP SLPANIELEC PLLGDVVICP EVLLAESAQF NKSLHAHWSL ILIHGVLHLL GYDHIKDEEA SIMQMLEAKL LAELGYANPY EVEENELE
Source	Mammalian cell
Target Names	ybeY
Protein Names	Recommended name: Probable rRNA maturation factor
Expression Region	1-158
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.