



Recombinant Legionella pneumophila subsp. pneumophila tRNA pseudouridine synthase A (truA)

Product Code	CSB-MP730557LDL
Abbreviation	truA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZVY6
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)
Purity	≥85% (SDS-PAGE)
Sequence	MRIALVVEYD GSQYHGWQAQ TGLHTIQQAV EFALSKVADS SISVVCAGRT DTGVHATNQV IHFDCEKDRS IRAWIHGANT FLPKDICVKW GKEMPENFHA RYSAVSRRYR YVIYNGAIRP GLLRGNVTWQ YRQLDHRLMQ QGGQCLLGEN DFTSFRSVEC QSNTPMRNIH QLQVTRHGDV VVLDITANAF LHHMVRNIAG VLIAVGSGKH PVGWVKDVLN AKDRKLGAEV APSYGLYLVQ VTYPKEFGLL QNNPGPLFLW EK
Source	Mammalian cell
Target Names	truA
Protein Names	Recommended name: tRNA pseudouridine synthase A EC= 5.4.99.12 Alternative name(s): tRNA pseudouridine(38-40) synthase tRNA pseudouridylate synthase I tRNA-uridine isomerase I
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.