



Recombinant Legionella pneumophila subsp. pneumophila Ribose-5-phosphate isomerase A (rpiA)

Product Code	CSB-EP730575LDL-B
Abbreviation	rpiA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZZB7
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)
Purity	≥85% (SDS-PAGE)
Sequence	MSELKIKAAK AAIAYIEDDM VIGVGTGSTV NFFIKELAAI KHKIEACVAS SKATEALLRA EGIPVIDLNS VQDLPIYVDG ADEVNERGEM IKGGGGALTR EKIVANVATQ FICIVDESKV VKRLGEFPVA VEVIPMARSF VARQIVKLGG DPEYREGFVT DNGNIILDVF NLSFSTPMAL EDSLNVIPGV VENGVFAKRL ADKVLVASAS GVNNLK
Source	E.coli
Target Names	rpiA
Protein Names	Recommended name: Ribose-5-phosphate isomerase A EC= 5.3.1.6 Alternative name(s): Phosphoriboisomerase A Short name= PRI
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.