



Recombinant *Methylococcus capsulatus* tRNA modification GTPase MnmE (mnmE)

Product Code	CSB-EP730606MFM-B
Abbreviation	mnmE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q602M5
Product Type	Recombinant Protein
Immunogen Species	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MAAPDRDTIT AIATPPGKGG VGIVRISGSN LGPVL DALLG RPPRPRYA EF RHFLDADGRA IDSGIALYFP APRSFTGENV LELHGHGGPV VLDLLLRRTL QLGCRLARPG EFSERAYLNG KDLAQAEAI ADLIDSSSTEE SARSAQRSLQ GEFSAHIHHL QECLVRLRTY VEA AIDFSDE DIDLLDDATL GHEITGLLDE LDTIDTKAHQ GALLREGLTT VIAGRPNVGK SLLNALAGR DLAI VTEIPG TTRDLLRESL QVGGLPLHIV DTAGLRDSED PIEREGIRRA RDALANADCI LLVCDARHTE AGDALPADLP ETIPLIRIFN KIDLTGAPAS LTVERETTVI HLSARTGEGV DLLRQEIRR AGYKKGTEGV FSARRRHLDA IRRARA AVAN ARLYLHTKTA ELLAEELRAA QKALGEITGE FTNEDLLNRI FSSFCLGK
Source	E.coli
Target Names	mnmE
Protein Names	Recommended name: tRNA modification GTPase MnmE EC= 3.6.-.-
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.