



Recombinant Mouse COUP transcription factor 1 (Nr2f1)

Product Code	CSB-MP730643MO
Abbreviation	Nr2f1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60632
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MAMVVSSWRD PQDDVAGGNP GGPNPAAQAA RGGGGGAGEQ QQAGSGAPHT PQTPGQPGAP ATPGTQGDKG QGPPGSGQSQ QHIECVVCGD KSSGKHYGQF TCEGCKSFFK RSVRRNLTYT CRANRNCPID QHHRNQCQYC RLKKCLKVGM RREAVQRGRM PPTQPMPGQY ALTNGDPLNG HCYLSGYISL LLRAEPTYTS RYGSQCMQPN NIMG IENICE LAARLLFSAV EWARNIPFFP DLQITDQVSL LRLTWSELFV LNAAQCSMPL HVAPLLAAAG LHASPMSADR VVAFMDHIRI FQEQVEKLKA LHVDSAEYSC LKAIVLFTSD ACGLSDDAHI ESLQEKSQCA LEEYVRSQYP NQPSRFGKLL LRLPSLRTVS SSVIEQLFFV RLVGKTPIET LIRDMLLSGS SFNWPYMSIQ CS
Source	Mammalian cell
Target Names	Nr2f1
Protein Names	Recommended name: COUP transcription factor 1 Short name= COUP-TF1 Alternative name(s): COUP transcription factor I Short name= COUP-TF I Nuclear receptor subfamily 2 group F member 1 V-erbA-related protein 3 Short name=
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.