



Recombinant *Methylococcus capsulatus* Enolase 1 (eno1)

Product Code	CSB-MP730652MFM
Abbreviation	eno1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q606T2
Product Type	Recombinant Protein
Immunogen Species	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	>85% (SDS-PAGE)
Sequence	MKARIESVQA MPIYDSRGVP TVRVRLTLD GTAATASIPS GASTGENEAV ELRDGDPKAH GGKGVLKAVS NVNDVIAPAV MGLEPWRQAE IDRLMIELD G TANKERLGAN AILGVSEAVA AAAAKSAGLP LYAYLGGAGQ ARLPIPMINI LNGGKHADSS LDFQEFMIMP VGAPSFAEAM RYATETFQAL KSILKAKGHA TSVGDEGGFA PQLQSNDEAC DLIVEAITKA GYEPGEDIAI ALDPAASSFH ENGVYRLTRS GQGDKTSSEM TDLYRRWIDK YPIVSIEDGL AENDWDGFRE HTAVLGEKIQ IVGDDLLVTN TRFIQRAIEE KSCNAVLIKL NQIGTVTETV EAIHLCKRAG WGFVISHRSG ETEDTFMADF AVAMGGGQIK TGSVCRSERM AKYNRLLEIE SDLGKAARFG R
Source	Mammalian cell
Target Names	eno1
Protein Names	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.