



Recombinant Mouse Immunity-related GTPase family M protein 1 (Irgm1)

| | |
|--------------------------|---|
| Product Code | CSB-BP730655MO |
| Abbreviation | Irgm1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q60766 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKPSHSSCEA APLLPNMAET HYAPLSSAFP FVTSYQTGSS RLPEVSRSTE RALREGKLLLE LVYGIKETVA TLSQIPVSIF VTGDSENGMS SFINALRVIG HDEDASAPTG VVRTTKRTE YSSSHFPNVV LWDLPGLGAT AQTVEDYVEE MKFSTCDLFI IIASEQFSSN HVKLSKIIQS MGKRFYIVWT KLDRDLSTSV LSEVRLQNI QENIRENLQK EKVYPPVFL VSSLDPLLYD FPKLRDTLHK DLSNIRCCPEP LKTLYGTYEK IVGDKVAVWK QRIANESLKN SLGVRDDDNM GECLKVYRLI FGVDESQQ VAQSMGTVM EYKDNMKSQN FYTLRREDWK LRLMTCAIVN AFFRLLRFLP CVCCCLRRLR HKRMLFLVAQ DTKNILEKIL RDSIFPPQI |
| Source | Baculovirus |
| Target Names | Irgm1 |
| Protein Names | Recommended name: Immunity-related GTPase family M protein 1 EC= 3.6.5.- Alternative name(s): Interferon-inducible GTPase 3 Interferon-inducible protein 1 LPS-stimulated RAW 264.7 macrophage protein 47 Short name= LRG-47 |
| Expression Region | 1-409 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.