



Recombinant Mouse Krueppel-like factor 4 (Klf4)

Product Code	CSB-MP730657MO
Abbreviation	Klf4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60793
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MRQPPGESDM AVSDALLPSF STFASGPAGR EKTLRPAGAP TNRWREELSH MKRLPPLPGR PYDLAATVAT DLESGGAGAA CSSNNPALLA RRETEEFNDL LDLDFILSNS LTHQESVAAT VTTSASASSS SSPASSGPAS APSTCSFSYP IRAGGDPGVA ASNTGGGLLY SRESAPPPTA PFNLADINDV SPSGGFVAEL LRPELDPVYI PPQQPPGG GLMGKVLKA SLTPGSEYS SPSVISVSKG SPDGSHPVVV APYSGGPPRM CPKIKQEAVP SCTVSRSLA HLSAGPQLSN GHRPNTHDFP LGRQLPTRTT PTLSPPELLN SRDCHPGLPL PPGFHPHPGP NYPPFLPDQM QSQVPSLHYQ ELMPPGSCLP EEPKPKRGRR SWPRKRTATH TCDYAGCGKT YTKSSHLKAH LRTHTGEKPY HCDWDGCGWK FARSDDELTRH YRKHTGHRPF QCQKCDRAFS RSDHLALHMK RHF
Source	Mammalian cell
Target Names	Klf4
Protein Names	Recommended name: Krueppel-like factor 4 Alternative name(s): Epithelial zinc finger protein EZF Gut-enriched krueppel-like factor
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.