



# Recombinant *Methylococcus capsulatus* Ribosomal RNA small subunit methyltransferase A (rsmA)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP730687MFM   |
| <b>Abbreviation</b>      | rsmA  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q60B77  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MRHVPRKRFG QNFLRDPGVI QEIVAAVGPA PSDRLVEIGP GEGVLTRELL<br>QSGACLEAIE LDRDLVAALK RRFAGVGRLR IHEGDAMKFD LRTIATGERL<br>RVVGNLPYNI STPLLHFLFD QIDVIEDMHF MLQKEVVDR L CAGAGDDHYG<br>RLSVMAALYC QAQHLFDVGP ECFHPQPKVV SAVVRLVPHA VPPDAGMVKQ<br>VSAVVAAFG QRRKTLRNAL KGLLDETAMV RAGIDPGARA EELSLADYVG<br>LSRQLNP   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | rsmA  |
| <b>Protein Names</b>     | Recommended name: Ribosomal RNA small subunit methyltransferase A EC=2.1.1.182 Alternative name(s): 16S rRNA (adenine(1518)-N(6)/adenine(1519)-N(6))-dimethyltransferase 16S rRNA dimethyladenosine transferase 16S rRNA dimethylase  |
| <b>Expression Region</b> | 1-257   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.