



Recombinant Mouse Myeloblastin (Prtn3)

Product Code	CSB-BP730722MO
Abbreviation	Prtn3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q61096
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	I VGGHEARPHS RPYVASLQLS RFPGSHFCGG TLIHPRFVLT AAHCLQDISW QLVTVVLGAH DLLSSEPEQQ KFTISQVFQN NYNPEENLND VLLLQLNRTA SLGKEVAVAS LPQQDQTLSQ GTQCLAMGWG RLGTQAPTPR VLQELNVTVV TFLCREHNVC TLVPRRAAGI CFGDSGGPLI CNGILHGVDV FVIRECASLQ FPDFFARVSM YVDWIQNVLR
Source	Baculovirus
Target Names	Prtn3
Protein Names	Recommended name: Myeloblastin EC= 3.4.21.76 Alternative name(s): Proteinase 3 Short name= PR-3
Expression Region	30-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>