



Recombinant Mouse SHC-transforming protein 3 (Shc3)

| | |
|--------------------------|---|
| Product Code | CSB-MP730724MO |
| Abbreviation | Shc3 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q61120 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSATRKS RAG DEPLRPPRG APHTSDQVLG PGVTYVVKYL GCIEVLRSMR SLDFSTR TQV TREAIRVCE AVPGAKGALK KRKPPSKMLS SILGKSNLQF AGMSISLTIS TASLNL RTPD SKQIIANHHM RSISFASGGD PDDTDYVAYV AKDPVNRRAC HILECCDGLA QDVIGSIGQA FELRFKQYLQ CPSKVPALQD RMQSLDEPWT EEGDGPDHP YNSVPTKMP PGGFLDARL KGRPHAPEAA QFAGKEQTY Y QGRHLGDTFG EDWQRAPTRQ GSLDIYSTAE GKTHMVPVGE TPTYVNTQPV PPQVWPAATS STESSPRKDL FDMKPFEDAL RNQPLGPMLS KAASVECISP VTPRAPDARM LEELNAEPWY QGEMSRKEAE ALLREDGDFL VRKSTTNPGS FVLTGMHNGQ AKHLLLVDPE GTIRTKDRVF DSISHLINYH LESSLPIVSA GSELCLQQPV ERKP |
| Source | Mammalian cell |
| Target Names | Shc3 |
| Protein Names | Recommended name: SHC-transforming protein 3 Alternative name(s): Neuronal Shc Short name= N-Shc SHC-transforming protein C Src homology 2 domain-containing-transforming protein C3 Short name= SH2 domain protein C3 |
| Expression Region | 1-474 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.