



# Recombinant Mouse Steroid hormone receptor ERR2 (Esrrb)

<b>Product Code</b>	CSB-MP730742MO
<b>Abbreviation</b>	Esrrb
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q61539
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSEDRHLGS SCGSFIKTEP SSPSSGIDAL SHHSPSGSSD ASGGFGIALS THANGLDSPP MFAGAGLGGN PCRKSYEDCT SGIMEDSAIK CEYMLNAIPK RLCLVCGDIA SGYHYGVASC EACKAFFKRT IQGNIEYNCP ATNECEITKR RRKSCQACRF MKCLKVGMLK EGVRLDRVRG GRQKYKRRLD SENSPYLNLP ISPPAKKPLT KIVSNLLGVE QDKLYAMPPN DIPEGDIKAL TTLCELADRE LVFLINWAKH IPGFPSLTLG DQMSELLQSAW MEILILGIVY RSLPYDDKLA YAEDYIMDEE HSRLVGLLDL YRAILQLVRR YKCLKVEKEE FMILKALALA NSDSMYIENL EAVQKLQDLL HEALQDYELS QRHEEPRRAG KLLLTLP LLR QTAAKAVQHF YSVKLQ GKVP MHKLFLEMLE AKV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Esrrb
<b>Protein Names</b>	Recommended name: Steroid hormone receptor ERR2 Alternative name(s): Estrogen receptor-like 2 Estrogen-related receptor beta Short name= ERR-beta Nuclear receptor subfamily 3 group B member 2
<b>Expression Region</b>	1-433
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a protein with similarity to the estrogen receptor. Its function is unknown; however, a similar protein in mouse plays an essential role in placental development.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.