



Recombinant Mouse Chromobox protein homolog 5 (Cbx5)

Product Code	CSB-EP730749MO
Abbreviation	Cbx5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q61686
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MGKKTkRTAD SSSSEDEEEY VVEKVLDRRM VKGQVEYLLK WKGfSEEhNT WEPEKNLDCP ELISEFMKKY KKMKEGEnNK PREKSEGNKR KSSFSNSADD IKSKKKREQS NDIARGFERG LEPEKIIGAT DSCGDLmFLM KWKDTDEADL VLAKEANVKC PQIVIAFYEE RLtWHAYPED AENKEKESAK S
Source	E.coli
Target Names	Cbx5
Protein Names	Recommended name: Chromobox protein homolog 5 Alternative name(s): Heterochromatin protein 1 homolog alpha Short name= HP1 alpha
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a highly conserved nonhistone protein, which is a member of the heterochromatin protein family. The protein is enriched in the heterochromatin and associated with centromeres. The protein has a single N-terminal chromodomain which can bind to histone proteins via methylated lysine residues, and a C-terminal chromo shadow-domain (CSD) which is responsible for the homodimerization and interaction with a number of chromatin-associated nonhistone proteins. The encoded product is involved in the formation of functional kinetochore through interaction with essential kinetochore proteins. The gene has a pseudogene located on chromosome 3. Multiple alternatively spliced variants, encoding the same protein, have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.