



Recombinant Mouse Serine/arginine-rich splicing factor 2 (Srsf2)

Product Code	CSB-MP730774MO
Abbreviation	Srsf2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62093
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SYGRPPPDV EGMTSLKVDN LTYRTSPDTL RRVFEKYGRV GDVYIPRDRY TKESRGFAFV RFHDKRDAED AMDAMDGAVL DGREL RVQMA RYGRPPDSHH SRRGPPPRRY GGGGYGRRSR SPRRRRRSRS RSRSRSRSRS RSRYSRSKSR SRTRSRSRST SKRSARRSK SKSSSVSRSR SRSRSRSRSR SPPPVSKRES KRSRSRSKSP KSPEEEGAVS S
Source	Mammalian cell
Target Names	Srsf2
Protein Names	Recommended name: Serine/arginine-rich splicing factor 2 Alternative name(s): Protein PR264 Putative myelin regulatory factor 1 Short name= MRF-1 Splicing component, 35 kDa Splicing factor SC35 Short name= SC-35 Splic
Expression Region	2-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.