



Recombinant Mouse Replication protein A 32 kDa subunit (Rpa2)

Product Code	CSB-EP730779MO-B
Abbreviation	Rpa2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62193
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MWNSGFESFT SSTYGGRRGGY TQSPGGFGSP TPSQAEKKS VRAQHVPCT ISQLLSATLT DEVFRIGDVE ISQVTIVGII RHAEKAPTNI VYKIDDMTAP PMDVRQWVDT DDASGENAVV PPETYVKVAG HLRSFQNKKS LVAFKIIPLE DMNEFTAHIL EVVNSHMMLS KPNSQASAGR PSMSNPGMSE SFNFGNNFM PANRLTVVQN QVLNLIKACP RPEGLNFQDL RSQLQHMPVP SIKQAVDFLC NEGHIYSTVD DDHFKSTDAE
Source	E.coli
Target Names	Rpa2
Protein Names	Recommended name: Replication protein A 32 kDa subunit Short name= RP-A p32 Alternative name(s): Replication factor A protein 2 Short name= RF-A protein 2 Replication protein A 34 kDa subunit Short name= RP-A p34
Expression Region	1-270
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.