



Recombinant Rat Cytochrome b-245 light chain (Cyba)

Product Code	CSB-EP730804RA-B
Abbreviation	Cyba
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62737
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	GQIEWAMWA NEQALASGLI LITGGIVATA GRFTQWYFGA YSIVAGVLIC LLEYPRGKRK KGSTMERCQK KYLTAVVKLF GPLTRNYYVR AVLHLLLSVP AGFLLATILG TVCLAIASVI YLLAAIRGEQ WTPIEPKPKE RPQVGGTIKQ PPTNPPPRPP AEVRKKPSEA EEEAASAGGP QVNPIPVTD E VV
Source	E.coli
Target Names	Cyba
Protein Names	Recommended name: Cytochrome b-245 light chain Alternative name(s): Cytochrome b(558) alpha chain Cytochrome b558 subunit alpha Neutrophil cytochrome b 22 kDa polypeptide Superoxide-generating NADPH oxidase light chain subunit p22
Expression Region	2-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Cytochrome b is comprised of a light chain (alpha) and a heavy chain (beta). This gene encodes the light, alpha subunit which has been proposed as a primary component of the microbicidal oxidase system of phagocytes. Mutations in this gene are associated with autosomal recessive chronic granulomatous disease (CGD), that is characterized by the failure of activated phagocytes to generate superoxide, which is important for the microbicidal activity of these cells.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.