



Recombinant Bacillus cereus Phosphate acyltransferase (plsX)

Product Code	CSB-EP730918BAAD-B
Abbreviation	plsX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q636H3
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MKIAIDAMGG DHAPKAVVLG AMKAIKEYSD LHITLVGKEE EIRQYLTSEE RITILHTDEK IESTDEPVRA VRRKKQASMV LAAQQVKDGV ADACISAGST GALMAAGLFV VGRMEGIERP ALSPTMPTVD GEGFVMLDVG ANVDAKPIHL YQYAVMGSVY AEKVRGIKNP RVGLLNVGTE DGKGNELSKQ VFTMLKDAPI NFVGNVESRD LLQGVADV VV CDGFTGNVAL KSLEGTALAL FSMLKEQLMS SFTSKLAAAV LKPKLMVLKD KMDYSEYGGA ALFGLKAPVI KAHGSSNDQS IFSAIRQTRE MVAKEVIPTI SSVMEKEPLQ
Source	E.coli
Target Names	plsX
Protein Names	Recommended name: Phosphate acyltransferase EC= 2.3.1.n2 Alternative name(s): Acyl-ACP phosphotransacylase Acyl-[acyl-carrier-protein]--phosphate acyltransferase Phosphate-acyl-ACP acyltransferase
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.