



Recombinant *Bacillus cereus* Glycerol kinase (glpK)

Product Code	CSB-YP730947BAAD
Abbreviation	glpK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63EX2
Product Type	Recombinant Protein
Immunogen Species	<i>Bacillus cereus</i> (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MKKYILSLDQ GTTSSRAILF NKKGEIVHSA QKEFTQHFPK PGWVEHNAQE IWGSILAVIA TCLSEADV KP EQIAGIGITN QRETAVVWDK TTGKPIYNAI VWQSRQTAEI CDELKEKGYG EMVREKTGLL IDAYFSGTKV KWILDNVEGA REKAENGDLL FGTIDTWLVW KLSGGKAHVT DYSNASRTL M FNIHDLQWDD ELLDM LTVPK SMLPEVRPSS EVYGETIDYH FFGQNVPIAG VAGDQQAALF GQACFGEGMA KNTYGTGCFM LMNTGEKAVA SEHGLLTIA WGIDGKVNYA LEGSIFVAGS AIQWLRDGMR MFKDASESEV YASRVESTDG VYVVPFVGL GTPYWDSEVR GAMFGVTRGT TKEHFIRATL ESLAYQTKDV LCAMEADSGI ELKTLRVDGG AVKNNFLMKF QSDILDVPVE RPVINETTAL GAAYLAGLAV GYWKNQDEIK EQWHMDKRFE PTMEAETSEE LYAGWKKAIE ATKAFK
Source	Yeast
Target Names	glpK
Protein Names	Recommended name: Glycerol kinase EC= 2.7.1.30 Alternative name(s): ATP:glycerol 3-phosphotransferase Glycerokinase Short name= GK
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.