



Recombinant *Xenopus laevis* STAM-binding protein-like (stambp)

Product Code	CSB-MP731019XBE
Abbreviation	stambp
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63ZM7
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MPEHSDASLP PEERIRALVL KGTSVEVNDD IPPKRYRSG VELIRMANVY SGEGSIENAF ILYNKYITLF IEKLPKHRDY KTANVPEKKE TLKKLKEIAF PKAEELKKEL HKRYKKEYEE YSEKQRKEEE ERARRLALQQ QLDAEKQRVA LLKQQQEQQE QVQAFEEMMR RKELEAERLR ILHQFSKDEP EAEPLGSPLI PGMNEPPVTP LLPSYGTVP HPPAVDRSLK PSSYGSNSSG VTSDGLRHVK IPRDVCKFL QLSNNTQRG VETCGILCGK LLQNEFTVTH VIVPKQSGGP DYCNTESEEE LFLIQDQQL ITLGIWHTHP TQTAFLSSVD LHTHCSYQMM LPESIAIVCS PKFQETGFFK LTDYGMKEIG ECRQKGFHPH CKEPPLFSAG GHVSVTEQDV TMMDLR
Source	Mammalian cell
Target Names	stambp
Protein Names	Recommended name: STAM-binding protein-like EC= 3.4.19.-
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.