



Recombinant Mouse Transcription factor E3 (Tfe3)

Product Code	CSB-BP731024MO
Abbreviation	Tfe3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64092
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MSHAAEPARD AVEASAEGPR AVFLLLEERR PAESAQLLSL NSLLPESGIV ADIELENILD PDSFYELKSQ PLFLRSSLP I SLQATPTTPA TLSASSSAGG SRTPAMSSSS SRVLLRQQLM RAQAQEQERR ERREQAAAAP FPSPAPASPA ISVIGVSAGG HTLSRPPPAQ VPREV LKVQT HLENPTRYHL QQARRQQVKQ YLSTTLGPKL ASQALTPPPG PSSAQPLPAP ETAHATGPTG SAPNSPMALL TIGSSSEKEI DDVIDEIIISL ESSYNDEMLS YLPGGTAGLQ LPSTLPVSGN LLDVYSSQGV ATPAITVSNS CPAELPNIKR EISETEAKAL LKERQKKDNH NLIERRRRFN INDRIKELGT LIPKSNPEM RWNKGTILKA SVDYIRKLQK EQQRSKDLES RQRSLEQANR SLQLRIQELE LQAQIHGLPV PPNPGLLSLT TSSVSDSLKP EQLDIEEEGR PSTTFHVSGG PAQNAPPQQP PAPPSDALLD LHFPSDHLGD LGDPFHGLGLE DILMEEEGMV GGLSGGALSP LRAASDPLLS SVSPAVSKAS SRRSSFSMEE ES
Source	Baculovirus
Target Names	Tfe3
Protein Names	Recommended name: Transcription factor E3
Expression Region	1-572
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.