



Recombinant Sulfolobus solfataricus Elongation factor 1-beta (ef1b)

Product Code	CSB-MP731038FPM
Abbreviation	ef1b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64214
Product Type	Recombinant Protein
Immunogen Species	Sulfolobus solfataricus (strain ATCC 35092 / DSM 1617 / JCM 11322 / P2)
Purity	≥85% (SDS-PAGE)
Sequence	TDVLVVLKV FPDSDEVNLD NLYTDISSKL PKEYKIIRKE TEPIAFGLNA LILYVQMPEQ TEGGTDNLEE VVNNIQGVSH AEVVGITRLG F
Source	Mammalian cell
Target Names	ef1b
Protein Names	Recommended name: Elongation factor 1-beta Short name= EF-1-beta Alternative name(s): aEF-1beta
Expression Region	2-91
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.