



Recombinant Danio rerio RNA binding protein fox-1 homolog 1 (rbfox1)

Product Code	CSB-MP731044DIL
Abbreviation	rbfox1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q642J5
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEEKGSKMVE QGNQDAPAPP ETMAQPFPSA QFAPPQNGIP AEYTPSHPHP TPDYSGQTPV AEHTLNLYTP AQTHSEPSGP DNSIQAVSGT ATQTDDSAQT DSQQQTQSSE ITEIKTQPKR LHSVNIPIFRF RDPDLRQMFQ QFGKILDVEI IFNERGSKGF GFVTFESSAD ADRAREKLHG TVVEGRKIEV NNATARVMTN KKTVPYANG WKLNPVVGAV YSPEFYAVPG FPYPAATAAA AAYRGAHLRG RGRTVYNTFR AAAPPHIPA YGGVVYQDGF YGADIYGGYT AYRYTQPATA TAAAYSDSYG RVYAADPYNH ALAPAATYSV GAMNAFAPLT DAKTRSHADD VGLVLSSLQA SIYRGGYSRF APY
Source	Mammalian cell
Target Names	rbfox1
Protein Names	Recommended name: RNA binding protein fox-1 homolog 1 Alternative name(s): Ataxin-2-binding protein 1 Fox-1 homolog A
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.