



Recombinant Mouse Forkhead box protein B1 (Foxb1)

Product Code	CSB-BP731079MO
Abbreviation	Foxb1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64732
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MPRPGRNTYS DQKPPYSYIS LTAMAIQSSP EKMLPLSEIY KFIMDRFPYY RENTQRWQNS LRHNLSFNDC FIKIPRRPDQ PGKGSFWALH PSCGDMFENG SFLRRRKRFK VLKSDHLAPS KPADAAQYLQ QQAKLRLSAL AASGTHLPQM PAAAYNLGGV AQP SGFKHPF AIENIIAREY KMPGGLAFSA MQPVPAAYPL PNQLTTMGSS LGTGWPHVYG SAGMIDSATP ISMTSGDYSY YGVPLKPLCH AAGQTLPAIP VPIKPTPAAV PALPALPAPI PTLLSNSPPS LSPTSSQTAT SQSSPATPSE TLTSPASALH SVAVH
Source	Baculovirus
Target Names	Foxb1
Protein Names	Recommended name: Forkhead box protein B1 Alternative name(s): Transcription factor FKH-5
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.