



Recombinant Bacteroides fragilis Holliday junction ATP-dependent DNA helicase RuvB (ruvB)

Product Code	CSB-MP731127BDQ
Abbreviation	ruvB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q650B4
Product Type	Recombinant Protein
Immunogen Species	Bacteroides fragilis (strain YCH46)
Purity	≥85% (SDS-PAGE)
Sequence	MEEDFNIRDH QLTSRERDFE NALRPLSFED FNGQDKVVDN LRIFVKAARL RAEALDHVLL HGPPGLGKTT LSNIANELG VGFKVTSGPV LDKPGDLAGV LTSLEPNDVL FIDEIHLSP VVEEYLYSAM EDYRIDIMID KGPSARSIQI DLNPFTLVGA TTRSGLLTAP LRARFGINLH LEYYDDDILS NIIRSAGIL DVPCSSQAAG EIASRSRGTP RIANALLRRV RDFAQVKGSG SIDTEIANYA LEALNIDKYG LDEIDNKILC TIIDKFKGGP VGLTTIATAL GEDAGTIEEV YEPFLIKEGF LKRTPRGREV TELAYKHLGR SLYNSQKTLF ND
Source	Mammalian cell
Target Names	ruvB
Protein Names	Recommended name: Holliday junction ATP-dependent DNA helicase RuvB EC= 3.6.4.12
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.