



Recombinant Orgyia pseudotsugata multicapsid polyhedrosis virus Uncharacterized 15.4 kDa protein (ORF73)

Product Code	CSB-YP731147OEU
Abbreviation	ORF73
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65366
Product Type	Recombinant Protein
Immunogen Species	Orgyia pseudotsugata multicapsid polyhedrosis virus (OpMNPV)
Purity	≥85% (SDS-PAGE)
Sequence	MWDTVRWQVL SSDEVEVAPD HRALAWRELI ANVARDTPLD YTFRTMLQRA DLENFDYNTP IVYNVKRKEL IVLNERMRAA LQRPVRRSDR TINANTAHVF LLFILAVLLT VLVAFSPWPD TRRVTAETRT I
Source	Yeast
Target Names	ORF73
Protein Names	Recommended name: Uncharacterized 15.4 kDa protein Alternative name(s): ORF73
Expression Region	1-131
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.