



Recombinant Protein grpE (grpE)

Product Code	CSB-EP731298BUI
Abbreviation	grpE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q661A2
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MEKKETKNDA EKNNKQDNKS TKSQKKNLN LVNSDKKITE LENEISNLKD LYLRKQAEFE NFRKRLEKEK DNFVKFANET IMKDVVNFLD NLERAINSSR KSKDFDNLIT GISMIENEIL SIFDKKYNLK KFGGENGFDF PSRHEAISIE EKEDLKNPEI VEVYQKGYCY NDRILRTAKV KVAQSKN
Source	E.coli
Target Names	grpE
Protein Names	Recommended name: Protein grpE Alternative name(s): HSP-70 cofactor
Expression Region	1-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.