



# Recombinant Elongation factor Tu (tuf)

<b>Product Code</b>	CSB-EP731304BUI-B
<b>Abbreviation</b>	tuf
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q661E5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Borrelia bavaricensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia bavaricensis)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAKEVFQRTK PHMNVGTIGH VDHGKTTLTA AISIYCSKLN KDAKALKYED IDNAPEEKAR GITINARHIE YETANRHYAH VDCPGHADIYI KNMITGAAQM DAAILLVAAD SGAEPQTKEH LLLAQRMGIK KIIVFLNKLD LADPELVELV EVEVLELVEK YGFSANTPII KGSAFGAMSN PEDPESTKCV KELLESMDNY FDLPERDIDK PFLLAVEDVF SISGRGTVAT GRIERGVIKV GQEVEIVGIK ETRKTTVTGV EMFQKILEQG QAGDNVGLLL RGVDDKDIER GQVLSAPGTI TPHKKFKASI YCLTKEEGGR HKPFFPGYRP QFFFRTTDVT GVVALEGKEM VMPGDNDIV VELISSIAMD KNVEFAVREG GRTVASGRIL EILE
<b>Source</b>	E.coli
<b>Target Names</b>	tuf
<b>Protein Names</b>	Recommended name: Elongation factor Tu Short name= EF-Tu
<b>Expression Region</b>	1-394
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>