



Recombinant Probable inorganic polyphosphate/ATP-NAD kinase (ppnK)

| | |
|--------------------------|---|
| Product Code | CSB-EP731310BUI |
| Abbreviation | ppnK |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q661V4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Borrelia bavariensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia bavariensis) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MKSKALLCIN TLKSGAGILG DDIKVYLETK HFVEVVLIDV SKPLLSFPRE NFLFLITLGG DGTVLLAVNL LLENKNVDIP IISINMGKVG FLADIKIEDF KKVIDRFFKN SLVINKKFL HVTVYKHGKD LISRYALNDI IIRSSLLNKM IHVDLKVNSE NFLSYKSDGI IVSTPTGSTG YSFSAGGPIL EAELEGFLT PISPHSVYNR SFVFSKLSKL SLSFSKEYFI APASIFLDGI NFGSFGVDVV FEFEISSQSL NFVSFCTDTF VKRLKNKLL |
| Source | E.coli |
| Target Names | nadK |
| Protein Names | Recommended name: Probable inorganic polyphosphate/ATP-NAD kinase Short name= Poly(P)/ATP NAD kinase EC= 2.7.1.23 |
| Expression Region | 1-279 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |