



Recombinant SsrA-binding protein (smpB)

Product Code	CSB-EP731318BUI
Abbreviation	smpB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q662X9
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBI) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MNTNTLLLEN KKAKFNYFIE EKVSCGIVLK GTEVKSIAK KLSFNNSFAS IKKEEFWLEN LHVSKYKEGN IFNHDELRRP KLLIKKQELQ RLKKFKEKEG YTLIPISFYL KKSIIKVEVG ICKGKKLYDK REILKQKSIK KDLSREIKYK
Source	E.coli
Target Names	smpB
Protein Names	Recommended name: SsrA-binding protein
Expression Region	1-150
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.