



Recombinant L-seryl-tRNA (Sec) selenium transferase (selA)

Product Code	CSB-MP731322YAH
Abbreviation	selA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q663U3
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MSAEPHPLYR QLPAIDRLN EPEMAPLLAE YGPVLLADTL RQLQAEAREY IGQFHTLADW CADWPAALRQ RLNQRQPALK PVFNLTGTVL HTNLGRAPLA ESAIAAVTDA MRSAVTLEYS LEGAGRHRD RAVADLLCAL TGAEDACIVN NNAAAVFLLL TVMAAGKQVV VSRGELVEIG GAFRIPDVMR QAGCDLVEVG TTNRTHLKDY RQAINENTGL LMKVHTSNYS IEGFTA AVSE QQLAALGQEC SIPTATDLGS GSLVDMTRYG LPAEPMPQQL IAAGVDLVTF SGDKLLGGPQ AGIILGKKQW IERLQQHPLK RALRADKMTL AALDATLRLY QQPDRLVEQL PSLRLLTRPA SEIAACAQRL LAPLIACYGT DFTLDIESCW SQIGSGSLPV DRLPSWALTF TPKDGRGSTL EALTARWRTL TKPVIGRVAD GRLWLDLRLCL EDEAALLREL AS
Source	Mammalian cell
Target Names	selA
Protein Names	Recommended name: L-seryl-tRNA(Sec) selenium transferase EC= 2.9.1.1 Alternative name(s): Selenocysteine synthase Short name= Sec synthase Selenocysteinyl-tRNA(Sec) synthase
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.