



# Recombinant Ribosomal RNA large subunit methyltransferase N (rlmN)

<b>Product Code</b>	CSB-MP731353YAH
<b>Abbreviation</b>	rlmN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q667Z6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Yersinia pseudotuberculosis serotype I (strain IP32953)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSEQLLTAST PIDAAPLSDN TVQTTAPATS KINLLDLNRQ QMREFFAEMG EKPFRADQVM KWMYHYCYDD FEQMTDINKG LRAKLQRVAE IRAPEVAEEQ RSV DGTIKWA IKVGDQQVET VYIPEADRAT LCVSSQVGCA LECKFCSTAQ QGFNRNLRVS EIIGQVWRAA KIIGSLKSTG TRPITNVMM GMGEPLLNLN NVVPAMDIMM DDFGFLSKR RVTLSTSGVV PALDKLGDMI DVALAISLHA PTDDIRDEIV PINRKYNIET FLAAVRRYLD KSKANGGRVT VEYVMLDHIN DSTEQAHLA ECLKDTPCKI NLIPWNPFPG APYGRSSNSR VDRFSKVLME YGFTTIVRKT RGDDIDAACG QLAGEVIDRT KRTLKKKMAG EPIAIKTV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rlmN
<b>Protein Names</b>	Recommended name: Ribosomal RNA large subunit methyltransferase N EC= 2.1.1.192 Alternative name(s): 23S rRNA (adenine(2503)-C(2))-methyltransferase 23S rRNA m2A2503 methyltransferase
<b>Expression Region</b>	1-398
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.