



# Recombinant Dihydroorotase (pyrC)

<b>Product Code</b>	CSB-EP731364YAH-B
<b>Abbreviation</b>	pyrC
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q669K2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Yersinia pseudotuberculosis serotype I (strain IP32953)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTAQPQTLKI RRPDDWHIHL RDDEMLSTVL PYTSEVFARA IVMPNLAQPI TTVASAIAYR ERILAAVPAG HKFTPLMTCY LTNSLDVKEL TTGFEQGVFT AAKLYPANAT TNSTHGVSDI PAIYPLFEQM QKIGMPLLIH GEVTDAAVDI FDREARFIDQ ILEPIRQKFP ELKIVFEHIT TKDAADYVLA GNRFLGATVT PQHLMFNRNH MLVGGIRPHL FCLPILKRST HQQALRAAVA SGSDRFFLGT DSAPHAKHRK ESSCGCAGVF NAPAALPAYA SVFEELNALQ HLEAFCALNG PRFYGLPVND DVVELVRTPF LQPEEIPLGN ESVIPFLAGQ TLNWSVKR
<b>Source</b>	E.coli
<b>Target Names</b>	pyrC
<b>Protein Names</b>	Recommended name: Dihydroorotase Short name= DHOase EC= 3.5.2.3
<b>Expression Region</b>	1-348
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.