



Recombinant Octanoyltransferase (lipB)

Product Code	CSB-YP731388YAH
Abbreviation	lipB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66DF4
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MMPRLQQHKI ILRQLGLQPY APVSQAMHNF TEFRTDTPD EIWLVEHQHV FTQQGAGKAE HVLMPGDIPV IQSDRGGQVT YHGPGQQVMY VMVDLKRAKI GVRQLVTAIE NTVIETLAHF NIDSHARPDA PGVYVEQQKI CSLGLRIRRG CSFHGLALNI AMDLEPFQRI NPCGYAGMQM TQVSALQPGV TVADVQPVLV REFTRQLGYP TAKLQPWSLS DYLLSSHSS SVL
Source	Yeast
Target Names	lipB
Protein Names	Recommended name: Octanoyltransferase EC= 2.3.1.181 Alternative name(s): Lipoate-protein ligase B Lipoyl/octanoyl transferase Octanoyl-[acyl-carrier- protein]-protein N-octanoyltransferase
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.