



Recombinant Acetylornithine deacetylase (argE)

Product Code	CSB-BP731417YAH
Abbreviation	argE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66G73
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MKMKLPPFIE LYRALIATPS ISAADSALDQ SNEALINLLA GWFADLGFRV EIQVPDPTRH KFNLLASIGE NENGEHGGGL LLAGHTDTPV YDEGRWTRDP FTLTEHDHKL YGLGTADMKG FFAFILDAVR DIDASKLTKP LYILATADEE TTMAGARYFA ANTQLRPDFA IIGEPTSLQP VRAHKGHISN AIRITGQSGH SSDPARGVNA IDLMHESITQ LMALRTTLQE RYHNPAFTIP YPTMNFHGHIN GGDAANRICA CCELHMDIRP LPGLTSLDLN ELMTEALEPV SQRWPGRLSI DELHPPIPGY ECPTDHHMVG VIEKLLGERT AVVNYCTEAP FIQQVCPTLV LGPGSINQAH QPDEFIDMAF IEPTRELIGQ LVDHFCQQK
Source	Baculovirus
Target Names	argE
Protein Names	Recommended name: Acetylornithine deacetylase Short name= AO Short name= Acetylornithinase EC= 3.5.1.16 Alternative name(s): N-acetylornithinase Short name= NAO
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.