



Recombinant Rat F-box only protein 39 (Fbxo39)

Product Code	CSB-EP731424RA-B
Abbreviation	Fbxo39
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66H10
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MDEDSEVTQP QDQSCWATLP DVCLRRVFWW LGDRDRSRAA LVC RKWNQIM YSADLWRYRT ITFSGRPSRV HASEFESALW YVKKFGRYLE HLEIKFLNPY NAVLTKKFQV TMRGLLSCLG KSNRLRLSLS IQHLELDRLV WRNSIRGSLI KSLSFLLKKM GKHLDHLSLK GARLTVEQGC HILNSLSYMR NENVASELNI EDDFSHHLAV YGSSQFNKAM ATFHNLTLT LNYNCISDEL LETLSENNAG TLRTMNIKCH VHDPHGQVWV GMSWAKLARQ ASNLKVNFFF ERVMKYERLA RILLQEIPVR SISLRSCYFS DPDWSMRPTL TDLLPTFRNT LQKLTFEFNN NHESLDEQLH LLILACRKLF YFKIWAFLDV KFVERILKSQ EEGQCSLRTL KVRIYTNRYE TNEEDRTLRE IYRKYRKLID SELNYFVIAY PMM
Source	E.coli
Target Names	Fbxo39
Protein Names	Recommended name: F-box only protein 39
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.