



# Recombinant Rat Protein FAM221B (Fam221b)

<b>Product Code</b>	CSB-MP731433RA
<b>Abbreviation</b>	Fam221b
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q66HD8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEADKTSEGP PATQDAKKHP SLTVPHAEDL PEATTPESPP HLLSSQEYFQ PPKVHSKTYP SLDLGTKDSE DLQDPSVSET PLKLDTANLP STPSQSSAQS VISLERAISK IPLDDFIYEN SISEVQEEPL PLSPTTDISK YEYPISDVQE EPLLLSPTTN ISKYEDSISD VQEEPLEDTP TADISETEYS ISDVQEEPLE PPPTADIPES EDYISDDHEE PLEPPPTADI PESSYEFISD KVPQTQVPES EPPFKHSVPE PSAQAKEAAS ADEEEAEEDD LTGGTSKTAAG SEHHARKK KEKRGCESRP VVPAKRAELV EMAMHRKQ FDDQVNDLFQ WEKNSSLKAI QTGIYIGWRC PHYLWDCFRI GDESKFCFGH LLKEHQIISD LSVPCSVSQC RCLMFCFIPS RPEEVGEFWL KKRATFDPKA WRAQCRCKHT HEEHAATGSH PCRHRGCYCN SFESNFLCAA CDRRWEEHET FFETEETKRR GKRPYGTNTA KNRHRPF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Fam221b
<b>Protein Names</b>	Recommended name: Protein FAM221B
<b>Expression Region</b>	1-507
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.