



# Recombinant Human Kazrin (KAZN)

<b>Product Code</b>	CSB-MP731476HU
<b>Abbreviation</b>	KAZN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q674X7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MMEDNKQLAL RIDGAVQSAS QEVTNLRAEL TATNRRLAEL SGGGGPGPGP GAAASASAAG DSAATNMENP QLGAQVLLRE EVSRLQEEVH LLRQMKEMLA KDLEESQGGK SSEVLSATEL RVQLAQKEQE LARAKEALQA MKADRKRLKG EKTDLVSQMQ QLYATLESRE EQLRDFIRNY EQHRKESEDA VKALAKEKDL LEREKWELRR QAKEATDHAT ALRSQLDLKD NRMKELEAEL AMAKQSLATL TKDVPKRHSL AMPGETVLNG NQEWVVQADL PLTAAIRQSQ QTLYHSHPPH PADRQAVRVS PCHSRQPSVI SDASAAEGDR SSTPSDINSP RHRTHSLCNG DSPGPVQKNL HNPVQSLED LEDQKRKKKK EKMGFGSISR VFARGKQRKS LDPGLFDDSD SQCSPTRQSL SLSEGEEQMD RLQQVELVRT TPMSHWKAGT VQAWLEVVMA MPMYVKACTE NVKSGKVLLS LSEDEDLQLGL GVCSSLHRRK LRLAIEDYRD AEAGRSLSKA AELDHHWVAK AWLNDIGLSQ YSQAFQNHV DGRMLNSLMK RDLEKHLNVS KKFHQVSILL GIELLYQVNF SREALQERRA RCETQNIDPV VWTNQRVLKW VRDIDLKEYA DNLTNSGVHG AVLVLEPTFN AEAMATALGI PSGKHILRRH LAEEMSAVFH PANSTGIREA ERFGTTPGRA SSVTRAGKEE NSSGLKYKAG RLPLGKIGRG FSSKDPDFHD DYGSLQNEDE GDDDPQSRLE QCRLEGYNLSL EVTNV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KAZN
<b>Protein Names</b>	Recommended name: Kazrin
<b>Expression Region</b>	1-775
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.