



Recombinant Energy-coupling factor transporter ATP-binding protein EcfA 1 (ecfA1)

Product Code	CSB-MP731507SRW
Abbreviation	ecfA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67JX4
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MPIVFEAVSH IYQPGTPFAW KALDDVSLTI PDGEFWGIIG PTGSGKSTLI QHMGLLRPT TGRVLVDGLD LSDRKTDLRK VRQRVGLVFQ YPEHQLFGET IFEDVAFGPR NMQLGAAEVE RRVMTALERV GLDPAMKDRS PFGLSGGQAR RVALAGVLAM EPRVLILDEP TAGLDPQGRR EILDLVGGFA RTGMTVVLVS HSMDDVAEYA DRVLVMHRGR VHMLGTPWEL FARRAELEAI GLGVPAVAL ADKLRERGS IPPDVVTMDE AVAAIRGVLA GRKGADRP
Source	Mammalian cell
Target Names	ecfA2
Protein Names	Recommended name: Energy-coupling factor transporter ATP-binding protein EcfA 1 Short name= ECF transporter A component EcfA 1 EC= 3.6.3.-
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.