



Recombinant Probable GTP-binding protein EngB (engB)

Product Code	CSB-EP731538SRW
Abbreviation	engB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67SJ6
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MGVHAEFVLS AVAQQQFPQD GLPEVALVGR SNVGKSSLIN ALVRNRKLAR TSNTPGRTQA LNFYRVWPQG KPRPEGEPQP DKDAGRTALS GPVLQAARES GAFYLVDMPG YGFARVSEAQ RREWALIEG YLLTRGALRG VLQIVDLRHP PTRDDVTMRE WIRHHRLPSL CVATKADKIG RTAWPRHRQV IARELGLDGD GEPLVLFSAE TGLGRDDVWR WIREHVQDWT FDM
Source	E.coli
Target Names	engB
Protein Names	Recommended name: Probable GTP-binding protein EngB
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.