



Recombinant UPF0042 nucleotide-binding protein STH186 (STH186)

Product Code	CSB-MP731540SRW
Abbreviation	STH186
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67T22
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MMDVDRSIRL VVTGMSGAG KTQALKYLED LGFFCVDNLP PSLMPKLAEL FGQTEGKVSRLALGIDIRGG RFFHEILGAL RQIAEIGVAY QILFMDASDE VLVRRYKETR RRHPLAAQGR VLDGIQRERR LLQELRGLAT FIIDTTHMTP ADLRKELNRR FGQDRESPHF HVNVVSFGFK HGAVLDADLV FDVRFLPNPH YVPDLQPLTG EDPAVVEYVM KWNVTQQFYR RLTGLIGFLL PHYVAEGKSL LTIAIGCTGG KHRSVCLANR LAHWIRERGY SVSVEHRDMP RPADRSDEEE QP
Source	Mammalian cell
Target Names	STH186
Protein Names	Recommended name: UPF0042 nucleotide-binding protein STH186
Expression Region	1-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.