



Recombinant Arabidopsis thaliana Probable aminotransferase TAT1 (At4g28420)

| | |
|--------------------------|---|
| Product Code | CSB-EP731549DOA |
| Abbreviation | At4g28420 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q67Y55 |
| Product Type | Recombinant Protein |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNHNSNLVLP SHQTETQTQD ETDISVWRFR GSDNAAKASS VTMRVIVYKL FDECSLDVKK PLLPLAHGDP SVYPCYRTSI LVENAVVDVL RSGKGN SYGP AAGILPARQA VADYVNRDLT NKVKPNDVFI TVGCNQGIEV VLQSLARPNA NILLPRPSYP HYEARAVYSG LEVRKFDLLP EKEWEIDLPG IEAMADENTV AMVIINPNNP CGNVYSYDHL KKVAETAKKL GIMVITDEVY CQTIFGDKPF VPMGEFSSIT PVITLGGISK GWIVPGWRIG WIALNDPRGI LKSTGMVQSI QQNLDITPDA TTIVQAALPE ILGKANKELF AKKNSMLKQN VELVCDRLKE IPCLVCNKKP ESCTYLLTKL KLPLLEDIED DMDFCMKLAK EENLVLLPGV ALGLKNWIRI TIGVEAQMLE DALERLNGFC KRHLKKTSS FQALSNGKI |
| Source | E.coli |
| Target Names | At4g28420 |
| Protein Names | Recommended name: Probable aminotransferase TAT1 EC= 2.6.1.- Alternative name(s): Tyrosine aminotransferase 1 |
| Expression Region | 1-449 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.