



Recombinant Arabidopsis thaliana Probable aminotransferase TAT1 (At4g28420)

Product Code	CSB-MP731549DOA
Abbreviation	At4g28420
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67Y55
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MNHNSNLVLP SHQTETQTQD ETDISVWRFR GSDNAAKASS VTMRVIVYKL FDECSLDVKK PLLPLAHGDP SVYPCYRTSI LVENAVVDVL RSGKGNISYGP AAGILPARQA VADYVNRDLT NKVKPNDVFI TVGCNQGIEV VLQSLARPNA NILLPRPSYP HYEARAVYSG LEVRKFDLLP EKEWEIDLPG IEAMADENTV AMVIINPNNP CGNVYSYDHL KKVAETAKKL GIMVITDEVY CQTIFGDKPF VPMGEFSSIT PVITLGGISK GWIVPGWRIG WIALNDPRGI LKSTGMVQSI QQNLDITPDA TTIVQAALPE ILGKANKELF AKKNSMLKQN VELVCDRLKE IPCLVCNKKP ESCTYLLTKL KLPLLEDIED DMDFCMKLAK EENLVLLPGV ALGLKNWIRI TIGVEAQMLE DALERLNGFC KRHLKKTSS FQALSNGKI
Source	Mammalian cell
Target Names	At4g28420
Protein Names	Recommended name: Probable aminotransferase TAT1 EC= 2.6.1.- Alternative name(s): Tyrosine aminotransferase 1
Expression Region	1-449
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.