



Recombinant Arabidopsis thaliana BTB/POZ domain-containing protein At1g55760 (At1g55760)

Product Code	CSB-YP731552DOA
Abbreviation	At1g55760
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q680K8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTDSAYRVDT ISRLAQWRIH NLSSTYRKS DPFKMGLWNW HLSVEKSKML LNVKLYPEVS NLTRENPVA SFALRVVST GERKALSHPE VIDKRIKTNE DFIWTIEVPL TGKIIDVEF LDLKVLQDS GELYSIWANG STENQSQVTA VTSLGRMLTE SIYTDITINA SDGSIGAHRA VLAARSPVFR SMFLHDLKEK ELSEINVLDM PLDACQAFLS YIYGNIQNED FLIHLALLQ AAEKYDIADL KEACHLSLLD DIDTKNVLER LQNAYLYQLP ELKASCMRYL VKFGKIFEIR DEFNIFMQCA DRDLISEIFH EVLSTWKGK
Source	Yeast
Target Names	At1g55760
Protein Names	Recommended name: BTB/POZ domain-containing protein At1g55760
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.