



Recombinant Arabidopsis thaliana Uncharacterized protein At1g03900 (At1g03900, At1g03901)

Product Code	CSB-MP731553DOA
Abbreviation	At1g03900, At1g03901
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q681Q7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSFEEEEEEE TFEHTLLVVR EVSVYKIPPR TTSGGYKCGE WLQSDKIWSG RLRVVSKDR CEIRLEDSNS GDLFAACFVD PGRRENSVEP SLDSSRYFVL RIDDGRGKYA FIGLGFAERN EAFDFNVALS DHEKYVRREK EKETGETSES DNHIDIHPAV NHRLKEGETI RINVKPKPTT NGTGMLSAAL SGTGKPKPLA LAPPPKAAGV TRSPLPPPPN DPVASRIASD GCKESRRNEP LSDLSQLKKN LPSTAGSGSS KSTGAASGWA AF
Source	Mammalian cell
Target Names	At1g03900
Protein Names	Recommended name: Uncharacterized protein At1g03900
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.