



Recombinant Human Olfactomedin-like protein 2A (OLFML2A)

Product Code	CSB-MP731560HU
Abbreviation	OLFML2A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68BL7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DSK VFGDLDQVRM TSEGSDCRCK CIMRPLSKDA CSRVRSGRAR VEDFYTVETV SSGTDCRCSC TAPPSSLNPC ENEWKMEKLLK KQAPPELLKLQ SMVDLLEGTL YSMDLMKVHA YVHKVASQMN TLEESIKANL SRENEVVKDS VRHLSEQLRH YENHSAIMLG IKKELSRLGL QLLQKDAAAA PATPATGTGS KAQDTARGKG KDISKYGSVQ KSFADRGLPK PPKEKLLQVE KLRKESGKGS FLQPTAKPRA LAQQQAVIRG FTYYKAGKQE VTEAVADNTL QGTSWLEQLP PKVEGRSNSA EPNSAEQDEA EPRSSERVDL ASGTPTSIPA TTTTATTTPT PTTSLLPTEP PSGPEVSSQG REASCEGTLR AVDPPVRHHS YGRHEGAWMK DPAARDDRIY VTNYYYGNL VEFRNLENFK QGRWSNMYKL PYNWIGTGHV VYQGAFYYNR AFTKNIKYD LRQRFVASWA LLPDVVYEDT TPWKWRGHSD IDFAVDESL WVIYPAVDDR DEAQPEVIVL SRLDPGDLSV HRETTWKTRL RRNSYGNCFV VCGILYAVDT YNQQEGQVAY AFDTHGTDA RPQLPFLNEH AYTTQIDYNP KERVLYAWDN GHQLTYTLHF VV
Source	Mammalian cell
Target Names	OLFML2A
Protein Names	Recommended name: Olfactomedin-like protein 2A Alternative name(s): Photomedin-1
Expression Region	28-652
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.