



Recombinant Human Metallo-beta-lactamase domain-containing protein 2 (MBLAC2)

Product Code	CSB-BP731562HU
Abbreviation	MBLAC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68D91
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSALEWYAHK SLGDGIFWIQ ERFYESGNRA NIWLVRGSEQ DVVIDTGLGL RSLPEYLYSS GLLQDREAKE DAARRPLLAV ATHVHFDHSG GLYQFDRVAV HHAEEAELAR GDNFETVTWL SDSEVVRTPS PGWRARQFRV QAVQPTLILQ DGDVINLGDR QLTVMHMPGH SRGSICLHDK DRKILFSGDV VYDGSLIDWL PYSRISDYVG TCERLIELVD RGLVEKVLPG HFNTFGAERL FRLASNYISK AGICHKVSTF AMRSLASLAL RVTNSRTSP
Source	Baculovirus
Target Names	MBLAC2
Protein Names	Recommended name: Metallo-beta-lactamase domain-containing protein 2 EC=3.-.-
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.